

# Fuji DRI-CHEM IMMUNO AU10V

## Referenzwerte



| Parameter | Einheit | Messbereich                                 | Hund*        | Katze*       |
|-----------|---------|---|--------------|--------------|
| T4        | µg/dL   | 0,5 - 8,0                                   | 1,3 - 2,9    | 0,6 - 3,9    |
|           | nmol/L  | 6,4 - 103,0                                 | 16,7 - 37,3  | 7,72 - 50,19 |
| TSH       | ng/ml   | 0,25 - 5                                    | < 0,5        | -            |
| Cortisol  | µg/dl   | 1,0 - 30,0<br>[Verdünnung: 1,0 - 50,0]      | 1,0 - 6,0 ** | -            |
|           | nmol/l  | 27,6 - 828,0<br>[Verdünnung: 27,6 - 1380,0] | 27 - 165 **  | -            |

### Quelle:

\* Herstellerangaben

\*\* ACTH Stimulationstest, Low/High Dose Dexamethason Test



scil animal care company GmbH  
Dina-Weissmann-Allee 6  
68519 Viernheim  
Tel.: +49 (0) 6204 78 90 222  
Fax: +49 (0) 6204 78 90 200  
Email: [info-de@scilvet.com](mailto:info-de@scilvet.com)  
[www.scilvet.de](http://www.scilvet.de)